



DentaComfort Pillow

Manufacturer:

AB Germa
Industrigatan 54-56
SE-291 36 Kristianstad
Sweden

Tel: +46 (0)44 12 30 30
Fax: +46 (0)44 10 31 79

E-mail: info@germa.se

www.germaprotec.com



Directions for use

Vacuum DentaComfort Pillow

The DentaComfort pillow equipment for dental practices is designed to help, support and stabilize the patients head and neck during treatment. It is designed to be used on most medical dentist chairs and offer ergonomic workplace during complex treatments and offer additional safety plus maximum of relaxing and comfort for the patient.

To create desire shape and support with DentaComfort pillow place the pillow under the patient´s head and neck and create the desired position, create vacuum inside the pillow by using a saliva sucker to make the pillow form stable in the desired position.

Storage

Store the pillow in room temperature and dry place.

Specification

DentaComfort pillow is made of PVC film, latex free outer cover with slide valve. Detached inner bag of Nylon filled with granulates of polystyrene. Inner bag is designed to control granulates and to prevent any loss of granulates if there arise a damage on the outer material.

Safe usage

1 Place the pillow under the patient´s head and neck and create the desired shape and position by lifting gently on the side with both hands.



2 Open the slide valve and suck out air with the saliva sucker until the desired stability is achieved by the vacuum. It is very quick to take out air.

Do not push too hard on valve towards the pillow, this to make sure you have free air flow.



3 Make sure to close the slide valve by sliding the valve to lock position and remove the saliva sucker. You have created the desired position; note hardness of the pillow is determent by how much air you take out from the pillow.



Warnings and precautions

Important, do not place the pillow in front of or on the face of a patient.

Make sure not to have the pillow close to any sharp tools, edges objects, strong heat or fire. Always avoid damaging the outside cover.

Inspect the pillow for any possible damages before use and make sure the pillow and valve is in good condition for use. Do not use the pillow if you have discovered any damages.

Cleaning and disinfection

Clean the pillow with water and soap after use. When called for disinfection it can be cleaned with most common disinfectants. Dry the pillow before next usage.

Operational room covers shall be washed according to tags and disinfections regulations.

Warranty

Warranty 1 (one) year from delivery date. No liability for outer violence or damages from improper handling.

Discard of product

DentaComfort pillows may be scraped with regular household garbage at the end of its duration. Material used in packing is environmentally friendly and can be recycled.

Accessories

Covers, operational room suitable cotton covers and single use cover (will be introduced)