## DIRECTIONS FOR USE

### VACUUM SPLINTS

**English:** Directions for use. **Germa Vacuum Splints – EasyFix PLUS**  
**Deutsch:** Bedienungsanleitung. **Germa Vakuumschienen – EasyFix PLUS**  
**Français:** Mode d'emploi. **Attelles à dépression Germa – EasyFix PLUS**  
**Español:** Instrucciones de uso. **Férulas de Vacío de Germa – EasyFix PLUS**  
**Svenska:** Bruksanvisning. **Germa Vakuumsplintar – EasyFix PLUS**  
**Italiano:** Istruzioni per l'uso. **Stecche a Depressione Germa – EasyFix PLUS**  
**Português:** Instruções de utilização. **Talas de Vácuo Germa – EasyFix PLUS**  
**Dansk:** Brugsanvisning. **Germa Vakuum Splints – EasyFix PLUS**  
**Nederlands:** Gebruiksaanwijzing. **Germa Vacuüm Spalken zijn – EasyFix PLUS**  
**Suomi:** Käyttöohje. **Germa Vacuum Splint – EasyFix PLUS**  
**Norsk:** Brukerveiledning. **Germa Vacuum Splint – EasyFix PLUS**

**Distributor:**

Ferno Norden AS  
Innlaget 298  
NO-3185 Skoppum  
Norway

Phone: +47 33 03 45 00  
Fax: +47 33 03 45 10

E-mail: [fno@fernonorden.com](mailto:fno@fernonorden.com)

**Manufacturer:**

AB Germa  
Industrigatan 54-56  
SE-291 36 Kristianstad  
Sweden

Tel: +46 (0)44 12 30 30  
Fax: +46 (0)44 10 31 79

E-mail: [info@germa.se](mailto:info@germa.se)

[www.germa.se](http://www.germa.se)

|             |                          |           |    |
|-------------|--------------------------|-----------|----|
| English:    | Directions for use       | Page      | 3  |
| Deutsch:    | Bedienungsanleitung      | Seite     | 6  |
| Français:   | Mode d'emploi            | Page      | 9  |
| Español:    | Instrucciones de uso     | Página    | 12 |
| Svenska:    | Bruksanvisning           | Sida      | 15 |
| Italiano:   | Istruzioni per l'uso     | Pagina    | 18 |
| Português   | Instruções de utilização | Página    | 21 |
| Dansk:      | Brugsanvisning           | Side      | 24 |
| Nederlands: | Gebruiksaanwijzing       | Bladzijde | 27 |
| Suomi:      | Käyttöohje               | Sivu      | 30 |
| Norsk:      | Brukerveiledning         | Side      | 33 |

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### 1. Warning and caution statements

#### CAUTION

- Check regularly that the vacuum pump works correctly and that the air can be evacuated from the Vacuum Splint so as to make it rigid.
- Avoid handling the Vacuum Splint close to sharp objects.
- Avoid handling the Vacuum Splint close to an open fire.
- Check regularly that the vacuum pump, the Vacuum Splint and the straps are not damaged or display signs of wear and tear. Repair any damages and change defective parts immediately.
- Check regularly that the Vacuum Splint has no cracks.

Recommendation: Usable life of 5 years

### 2. Intended use

The **EasyFix PLUS Vacuum Splints** are intended to stabilise injured patients during transport. They are suited for fixation of patients with hand, arm, leg, shoulder and hip injuries. The Vacuum Splints are available as full leg-, lower leg-, arm- splints.

### 3. Specifications

The EasyFix Vacuum Splints are in conformity with ISO 6940:1984

The EasyFix Vacuum Splints are in conformity with Council Directive 93/42/EEC concerning Medical Device

**Temperature range:** -30 °C to +70 °C

#### Parts/Materials

Cover material: PVC coated polyester fabric (flame retardant)

Inner bag: Polyester

Granules: Polystyrene granules (flame retardant)

| Dimensions                          | Length (cm) | Width (cm) | Height (cm) | Weight (kg) | Capacity for hand pump   |
|-------------------------------------|-------------|------------|-------------|-------------|--------------------------|
| <b>EasyFix PLUS Full Leg Splint</b> | 135         | 72         | 3           | 1,55        | -450hPa at c. 40 strokes |
| <b>EasyFix PLUS Leg Splint</b>      | 97          | 62         | 3           | 1,15        | -450hPa at c. 25 strokes |
| <b>EasyFix PLUS Arm Splint</b>      | 74          | 39         | 3           | 0,75        | -450hPa at c. 14 strokes |
| <b>EasyFix PLUS Child Splint</b>    | 40          | 30/36      | 3           | 0,45        | -450hPa at c. 9 strokes  |

### 4. Description

The Vacuum Splint consists of an outer cover made of hard-wearing PVC coated fabric with an innerbag of polyester filled with polystyrene granule. The inner bag is sewn in channels enabling the granule to spread evenly.

When the air is evacuated, the granule is moulded into a tight substance making the Vacuum Splint rigid and unable to shape. Before and while the air is evacuated the Vacuum Splint may be moulded according to the shape of the patient's limb in order to obtain optimum support of the patient.

The Vacuum Splint is equipped with straps with Velcro fastenings and for Full Leg splint - 2 extension straps. A head strap is optionally available for the Leg Splint and the Full Leg Splint.

## 5. Operating instructions

Full Leg Splint – see 1, Leg Splint – see 2, Arm Splint – see 3, Forearm Splint – see 4, Wrist Splint – see 5.  
(Pictures page 36)

1. Put the splint on an even surface near the patient.
2. Mould the splint so that it is even or, if necessary, in an appropriate state for the emergency.
3. Close the valve by turning it to the right (clockwise).
4. Place the leg/lower leg/arm/forearm/wrist in the medically desired position on the splint.
5. Mould the splint around the leg/arm so that the fracture becomes stable and the unit is comfortable for the patient.
6. Attach the Velcro fasteners and tighten them.
7. Lift the foot part of the splint to a supporting angle of the foot and fasten the strap crosswise to the leg.
8. Connect the pump. Evacuate the air, and continue to mould the splint. Tighten the Velcro fastenings once more.
9. The patient is ready to be transported.

For patients with short legs the foot part should be folded.

By using the extension straps it is possible to support fractures in hip or shoulder.

In special cases, the Full Leg Splint can be used as a vacuum mattress for a small child. For this you need the optional head strap.

In special cases, the Leg Splint can be used as a vacuum mattress for a baby. For this, you need the optional head strap.

## 6. Cleaning and disinfecting

After each use the Vacuum Splint should be wiped with a moist cloth and dried. If dirty the splint may be washed with soapy water and rinsed with clean water.

Allow the Vacuum Splint to dry completely before storing it.

If the Vacuum Splint needs disinfection first clean it carefully and let it dry. Wipe the complete splint with a solution containing 70% alcohol or use the disinfectant called LYSETOL AF.

Test the Vacuum Splint at regular intervals to make sure that it has the necessary tightness.

If leaks are found; see section 7 for repair instructions.

## 7. Maintenance and repair

At a regular basis the Vacuum Splint and the pump must be checked for visible damages.

At the same time the splint must be tested as to tightness. Furthermore the suction force of the vacuum pump must be checked according to Section 8.

Leaks can be found on the Vacuum Splint by visual inspection or by filling the splint with air covering it with soapy water all over the surface and the welds and watching out for air bubbles from the splint.

If leaks have been found on the splint, mark the leak, let any air out, clean the area around the leakage and let it dry.

Find a patch in a suitable size and put glue on the patch as well as on the leakage on the splint. Wait for the glue to become dust dry (approx. 2 minutes) and place the patch on the leakage. Firmly press the whole area of the patch establishing good contact and press out all air under the patch.

The repair will have full strength after 24 hours.

## 8. Test of the Vacuum Splint and Pump

Straighten out the Vacuum Splint, **close the valve** by turning it clockwise and connect the pump.

Evacuate the air until the Vacuum Splint feels rigid and leave it for about 2 hours; if the splint is as rigid as it was when you started it may be considered tight. Open the valve by turning it counter-clockwise.

If the Vacuum Splint is not airtight, see section 7 for repair or check the valve for leakage. Try to insert a new valve and repeat the test.

If the pump is not emptying the Vacuum Splint check the pump by closing the connection with your fingertip and pump a few times.

A negative pressure should be created in the tubing that can be felt at your fingertip. If the pump does not function do not repair it but return it for service to your dealer or to an Germa service centre.

## 9. Storage

Loosen the valve by turning it counter-clockwise and put the Vacuum Splint on an even surface.

When the granules have been spread evenly, fold the Vacuum Splint in its initial position.

Repack the Vacuum Splint so that it fits into your storage area. The Vacuum Splint is often positioned folded in a bag.

## 10. Ordering information – EasyFix Splint

|             |                               |
|-------------|-------------------------------|
| 272 111 000 | EasyFix PLUS Full Leg Splint  |
| 272 112 000 | EasyFix PLUS Leg Splint       |
| 272 113 000 | EasyFix PLUS Arm Splint       |
| 272 118 000 | Child EasyFix PLUS Arm Splint |

## 11. Spare parts

|             |                           |
|-------------|---------------------------|
| 15000000112 | Repair kit – EasyFix PLUS |
| 12633800000 | Multifit ValvS            |

## 12. Accessories

|             |   |
|-------------|---|
| 272 000 701 | Germa Vacuum hand vacuum pump   |
| 271 000 701 | Germa Vacuum foot vacuum pump   |
| 15271000710 | Carrying bag for splints  |
| 272 114 000 | Set of splints – EasyFix PLUS I (Arm Splint, Leg Splint, repair kit, hand pump and bag)                   |
| 272 115 000 | Set of splints – EasyFix PLUS II (Arm Splint, Leg Splint, Full Leg Splint, repair kit, hand pump and bag) |

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








**Report no. :**

(To be assigned in by  
Manufacturer)

**Please print or type or fill in by MS Word, all sections a. – i. Fax or e-mail to your supplier. Please use repair order when reporting repairs not associated with complaints.**

For claims please specify :  
Purchase Date

**Description of product**

**Description of Incident/Problem:** 1) What, 2) How?

When? In connection with: 1) transport, 2) unpacking, 3) preparing, 4) test, 5) cleaning, 6) normal use, 7) training, 8) disassembly, 9) assembly, 10) storage, 11) other:

**Info on the product, prior to problem,**

Storage: 1) stock, 2) site of use, 3) packed, 4) unpacked, 5) softpack, 6) other:

Used cleaning and disinfecting agents/Sterilizing method:

Pattern of use: 1) stand-by for emergency, 2) routine use, 3) stand by, 4) other:

Environment: 1) Indoor, 2) Outdoor, 3) Heat, 4) Cold, 5) Humidity, 6) Sunlight, 7) Gases, 8) other:

---

Which **actions** have been: a) taken, b) planned by the reporter or user at the site of incident?

1) examination of the product on site, 2) replacement of the product, 3) other:

---

Suggested actions by the distributor/complainant

---

---

**Returning of samples/faulty product(s),**

- 1) Already shipped/Will be shipped to \_\_\_\_\_,
- 2) Pending request from manufacturer,
- 3) Product/parts not available:

---

Date:

---

**We acknowledge the receipt of this report by manufacturer and will start processing it immediately. Please use the assigned complaint report number for further reference.**

**Your contact person in this matter is:** \_\_\_\_\_

Name:

Direct Tel.

Direct Fax.

Init/Date:

---