

## Step by Step Evacuation Set

Step by Step is quick and easy to apply to the mattress and includes a restraint system for the patient to be secured on the mattress.

This solution deliver a secure evacuation of laying patients with ergonomics and less hassle for the rescue teams.

### Product information

- Uses little space
- Weighs only 1,25 kg
- Can be used on all our transfer mattresses
- Can be used on weights up to 250 kg
- Easy to understand user manual
- Reinforced straps with adjustable handles for head and foot end
- Two patientbands for fixation
- Sewn bag prevents the patient to slip out

### Productno.

GE10001

Step by step evacuation set

