



## EC Declaration of conformity no.: 111110-002

|                   |  |
|-------------------|--|
| Product Name:     | EasyFix Vacuum Mattress Ready2go   |
| Germa Article No: | 23271501001  |
| Description:      | <p>Vacuum mattress, green/green, 208 x 130/72 x 4,5 cm, 10 handles (5 on each side) and valve at head-end. It has a pillow construction and is adjusted to fit restraint systems.</p> <p>Vacuum Mattress is a tool for fixation, immobilisation and transport of patients with suspected back- neck- or head injuries. The product consists of an airtight shell with handles and fixation straps. Innerbag in polyester, polystyrene granules and a vacuum valve.</p> <p>The mattress corresponds to European standard EN 1865 and EN 1789. Note that outside the ambulance the patient shall be secured by the belt system of the vacuum mattress. Inside the ambulance the patient shall be secured by the restraint system of the stretcher.</p> <p>Ready2go includes:<br/>EasyFix Vacuum Mattress, Vacuum foot pump, Repair kit EasyFix, Bag for mattress and Belt set buckle EasyFix</p> |
| Accessory:        |  |

Identification: All products with serial numbers issued after LOT number: **030704**, date: **2003-07-04** (yyyy-mm-dd).

We the manufacturer hereby declare that this product is in conformity with the requirements in:

- Annex I and Annex VII and Annex IX in Council directive 93/42/EEC. Equipment class: I, non-sterile, no measuring function.
- Swedish Medical Devices Act 1993:584

Signed in Kristianstad, Sweden on behalf of AB Germa by:

Position: Quality Manager  
Name: Johan Attelid

Date: 10/11-2011

Sign: 