



Easy Transfer Comfort Mattress



Easy Transfer Mattress is designed in co-operation with professional users focusing on ergonomics patient handling, patient safety and comfort for patient.

Easy Transfer Mattress has been tested by trauma teams to secure efficient handling plus quick, safe and easy transfer of patients including heavy duty patients.

Easy Transfer Mattress approved according to EN 1865 and EN 1789

Art no: 271 030 000

Features and Benefits:

Safety

Easy to adjust safety belts according to standard EN 1865 and EN 1789

Comfort

Extra comfort by using mix of Sensus memory foam 20 mm and comfort 40 mm V-foam which support patients body weight and deliver a safe and comfort feeling.

Ergonomics

Strong handles to support lift, carry or slide of patient on the mattress using 10 strong handles. It is possible to lower the head for upper airway management. Carry capacity 250 kg.

Designed to fit on top of Mondial stretchers and other most common stretcher designs.

- **Hygienic** full welded easy to clean with most common disinfectants
- **Sliding material / reducing friction** material under the mattress, easy to slide patient over to e.g. hospital bed.
- *Fits together with Germa Step by Step carry Kit



Accessories: Step by Step
Art no: 15 000 000 117
15 000 000 023

Specification - Easy Transfer Comfort Mattress

PVC- Polyester fibre flame retardant
PP-plate and grass fibre material to reduce friction
Foam Sensus — V-shaped (material that deliver super comfort and side support to the

Dimensions and weight

Length: 185 x 50 x 6 cm
Weight: 4,7 kg
Handles: 10

Ordering Specification

Easy Transfer Comfort Mattress

Art. No.

271030000

