

Germa Pillow for Stretcher/Transfer Mattress

The Germa Pillow is intended for use together with the Germa Stretcher/Transfer Mattress. The flexible pillow that is filled with granulate allows the patient to feel comfortable. Please note that it is not a vacuum product. The Pillow is attached to the stretcher itself with two Velcro straps.

Fits optional stretchers.

Can be used with optional Stretcher/Transfer Mattresses.

Material

Outer material: Polyester fibre, PVC coated (flame retardant)

Granules: Styroporgranules (flame retardant)



Product information

Productno.

Accessories



GEAS160

Germa Vacuum Hand
Pump