

## Mondial Flat Mattress

Mondial Mattress is developed in co-operation with Ferno International focusing on safety and comfort for the patient. This mattress is designed with help of trauma teams to secure efficient handling with safety and comfort for the patient. Mattress is designed to be compliant to the European standard EN 1865 and EN 1789.

**Materials:**

Padova PVC- Polyester fibre flame retardant, compliant to EN 1865, polyether.



### Product information

- Hygienic, all welded product, easy to clean with most common disinfectants
- Side support via V-shaped foam in different density in the mattress to create safe feeling and comfort.
- Designed to fit on top of Ferno Mondial stretchers, fixation with plastic hook. Important to use stretcher restraint system according to European standard EN 1789
- Velcro, hooks and loops adhesive opportunity for head part fixation of the mattress on to the stretcher.

**Productno.**

F2801

Mondial Flat Mattress

