

Transfer Mattress Saver / Blue line

Transfer Mattress is intended to be used as very comfortable stretcher mattress, its flexibility ensuring good fit on the stretcher. Furthermore, it can facilitate patient transport and/or transfer from one stretcher to an hospital bed or to another stretcher. The Transfer Mattress can be used as a soft stretcher and be carried with you to narrow environments to support lifting and carry of a patient using the handles, carry capacity of 250 kg.

The six handles on each side can be folded in under the mattress when laying on top off a stretcher and be folded out when needed.

Material:

Outer material: Polyester fibre, PVC coated (flame retardant)

Mattress alternatives:

Stretcher/Transfer Mattress 12 (6 handles on each side)

Stretcher/Transfer Mattress 16 (6 handles on each side and additional 2 handles at the ends)



Product information

- Durable sliding surface on the reverse of the mattress, glassfibre fabric
- Allows easy sliding of the mattress with patient on the top e.g. from stretcher/cot over to hospital bed.
- An evacuation tool used together with Step by Step unit.
- Six handles that allow easy lifting of patient up to 250 kg(500 lbs.)
- Approved according to EN1865

Productno.

271708000

Saver Stretcher/Transfer Mattress 12



Accessories



GE10001

Step by Step Evacuation
Set