



Safe Transfer II Mattress

Safe Transfer Mattress is developed in co-operation with our customers focusing on hygienic, ergonomic patient handling, patient safety and comfort.

Safe Transfer II Mattress has been tested together with medical trained trauma teams to secure efficient handling plus quick and easy transfer of patients also heavy duty patients.

Safe Transfer II Mattress has been tested and approved according to EN 1865 and EN 1789

Art no: 271 029 001

Features and Benefits:

Safety

Easy to adjust safety belts according to standard EN 1865 and EN 1789

Comfort


Extra comfort by using Sensus memory foam which support patients body weight and deliver a safe and comfort feeling.

Ergonomics

Easy to lift, carry or slide over patient just on the mattress with 12 strong handles. It is possible to lower the head for upper airway management. Carry capacity 250 kg.

Designed to fit on top of Allfa stretchers.

- Hygienic full welded easy to clean with most common disinfectants
- Sliding material under the mattress, easy to slide patient over to e.g. hospital bed.
- *Fits together with Germa Step by Step carry Kit



Specification - Safe Transfer II Mattress

PVC and PU - Polyester fibre flame retardant
4,5 cm polyether + 1,7 cm Sensus (material that forms easy and deliver super comfort)

Dimensions and weight

Length: 190 x 50 x 6 cm
Weight: 5,3 kg
Handles: 12

Ordering Specification

Safe Transfer II Mattress

Ferno—Germa label

Art. No.

21271029001

