



AB GERMA®

Carry / Slide Sheet

Heavy Duty Patient Handling.

Application

Ferno Carry Sheet Heavy Duty contains a PVC fabric material welded with 20 carry handles on all sides of the sheet.

Heavy Duty Carry Slide sheet designed to support and help the paramedics to transfer, slide carry a patient from narrow environments and/or bed or floor on to an ambulance stretcher and in next step from the stretcher to a hospital bed. The purpose is to support handling of heavy-duty patient with safety and improved ergonomic, appreciated ergonomic tool for heavy-duty patient handling.

Features and Benefits

- 20 handles on each side make it easy to find balance for lifting and carry a patient on board
- Sliding surface on the back of the sheet, reduction of friction when sliding
- Easy to clean fully welded product
- Durable and strong material and design with a carry capacity of 380 kg
- Slide sheet easy to store in the ambulance, accessible when needed.



Specifications

Approved according to European Standard
EN 1865
CE—marked
X-ray and MRI translucent

Length: 205 cm
Width: 143 cm
Weight: 4,3 kg

Material:

PVC—polyester fabric material (flame retardant)
Rope: Nylon
Fiberglass reverse

Ferno Carry / Slide sheet Heavy Duty

Art no: 005 000 005

**205 cm long, 143cm wide
Weight 4,3 kg**