



**AB GERMA®**

## EasyLift Shoulder – carry support

Easy Lift carry support EasyLift is a tool which is designed to ease the burden especially the hands and arms for heavier lifting.

The advantage is that with rapidly occurring situations it is possible to let go of the handles as with traditional carrying.

### Product information

- Easy to use, easy to understand
- Distributes the load
- Immediate disengagement at the unexpected event
- Easy to clean
- The friction lock support you
- Can be worn under a jacket / vest, making it immediately available
- Designed and Developed by ambulance personnel with extensive experience.
- Easy Lift is made of reinforced fabric with plastic / metal runners and a friction PVC part.
- Easy Lift is set individually

**Productno:** 15 000 000 006 EasyLift Shoulder green  
15 000 000 016 EasyLift Shoulder black

