

EasyFix PLUS Vacuum Splints

Application

Germa EasyFix PLUS line of splints:

The EasyFix PLUS is built on the original high quality, long-life model.

Easy Fix PLUS Vacuum splints are intended and designed to deliver premium stabilisation to injured patients during transport. They are suited for fixation of patients with hand, wrist, arm, leg, shoulder and/or hip injuries. The splint moulds around the individual fracture giving an absolute optimal and very comfortable support to the patient.

Designed and constructed to reduce needed preparation and is very reliable and a product you can trust, world class equipment.

- Anti-Bacterial hooks and loops
- Improved hygienic— BioSafe straps
- Flexible to positioning of the splint— premium immobilisation

Material

Outer material: Polyester fibre, PVC coated (flame retardant)

Inner bag: Polyamide

Granules: Styropor granules (flame retardant)



Arm Splint



Leg Splint



Full Leg Splint

Features	EasyFix
Unique inner bag with channel system, ensuring granulates are evenly dispersed. <i>Quick and easy to use.</i>	X
Robust weldings. <i>Trustful</i>	X
PVC-coated fabric that is easy to clean and disinfect.	X
Full Leg splints with extra Velcro for strap allowing the splint to be used as vacuum mattress for babies and children.	X
Full leg splint is provided with additional Velcro straps to ensure one size fits all.	X
Velcro is sewn onto a piece of PVC coated fabric to reinforce the strength of the Velcro straps.	X



Ordering Specification			Art. No EasyFix
Full Leg Splint	Immobilisation of full leg, hip or shoulder injuries to adults. Also works as a vacuum mattress for small children or babies.	136 x 73 x 3 cm Weight: 1.55 kg	272 111 000
Leg Splint	Immobilisation of lower leg on adult and full leg, hip or shoulder on child. Also works as a vacuum mattress for babies.	98 x 61 x 3 cm Weight: 1.15 kg	272 112 000
Arm Splint	Immobilisation of arm, forearm and wrist fractures.	74.5 x 39 x 3 cm Weight: 0.75 kg	272 113 000
Splint Set I	Arm, Leg, Repair Kit, Pump and Bag.	75 x 37 x 18 cm Weight: 3.10 kg	272 114 000
Splint Set II	Arm, Leg, Full leg, Repair Kit, Pump and Bag.	75 x 37 x 18 cm Weight: 4.65 kg	272 115 000
Vacuum Pump	One-hand-operated Vacuum Pump.	Weight: 0,6	272 000 701