



AB GERMA®

SPIDERSTRAP BACKBOARD X



PRODUCT SHEET

X-Spider is a outstanding product for fixation of a patient on top of a Scoop stretcher or Spine Board. It is developed in cooperation with professional Swedish ambulance personnel.

X-Spider is developed to improve immobilisation of the patient by using x shaped fastening straps over the sternum and pelvis, the X connection is placed on the sternum so it does not affect the breath of the patient. The X shape also makes it easier for the ambulance personel to adjust the straps to accomodate the size and length of the patient.

X-Spider lock with Velcro alternative plastic buckle and supplied in a handy storage bag

SPIDERSTRAP BACKBOARD X



PRODUCT SPECIFICATIONS

PRODUCT INFO

- 4 diagonal fastening straps and one central X shaped fastening strap.

Warranty period: 2 years

MATERIAL

- Polyester fibre
- Velcro or plastic buckles

SPECIFICATIONS

- Length of straps:
 - X-strap Chest 175 cm
 - Pelvic 162 cm
 - Upper leg 146 cm
 - Lower leg/foot 123 cm
- Weight 690 g
- Loading capacity: 380 kg

ORDER INFORMATION

Art no: 271000734

Spiderstrap backboard Velcro X

Art no: 271000735

Spiderstrap backboard Buckle X

Art no: 271000736

Color Coded Spiderstrap backboard Velcro X

BENEFITS

- Twelve-point, six strap system used to secure patient to backboards
- Size is adjustable from pediatric to adult Premium immobilisation of the patient.
- Limited pressure on the torso- easier for patient to breath.
- Option of X fixation of hips/Pelvis.
- Easy to release the patient from the scoop board/stretcher when arriving at the hospital.
- Can be ordered with a colored coded version for a more intuitive use.

Visiting address
AB Germa
Industrig. 54-56
SE-291 36 Kristianstad, Sweden

Post address
AB Germa
Box 2127
SE-291 02 Kristianstad, Sweden

www.germa.se
Phone: +46 (0)44 12 30 30
info@germa.se