

## Child EasyFix vacuum mattress

Germa offer new design Child EasyFix Vacuum mattress.

Child EasyFix is developed from our earlier experiences from the original high quality materials traditionally used and producing Germa world class vacuum mattresses.

Child EasyFix vacuum mattress using body shape design to optimize the vacuum mattress capabilities and limited storage volume. Intended immobilisation use for a child within the length of 90 cm up to 160cm this translate to a child between the age of approx. 5 years to 14 years. Child EasyFix unique construction which deliver efficient fixation of upper back, neck and head. Reinforced with base plate that cover cervical spine part. Together with fixation restraint system using X-system to support and create most effective immobilisation.

Recommended patient length capacity 90 up to 160 cm.

Note! usable outside this range if effective immobilisation is deemed as having been achieved

### Specifications

Maximum shrinkage on length: 1%

Maximum shrinkage on width: 3%

Approved according to European Standard

EN 1865 and EN 1789

### Material:

Outer material:

Polyester fibre, PVC-coated (flame retardant)

Inner bag: Polyester

Granules: Styropor granules (flame retardant)

Wood board upper part: Birch plywood

### Weight:

4,3 kg (105 lb) without fixation straps

4,7 kg with fixation straps

5,8 kg total weight including fixation straps + footrest pump + bag



### Product information

- Head, neck and upper back part improved with wood part to secure stronger immobilisation
- Improved chambers system
- Head strap system
- Compliant to EN 1789 and EN 1865
- Carry capacity of 150 kg, 6 handles
- Light weighted and small storage

### Productno.

271 504 000

Child EasyFix Vacuum Mattress 5 handles on each side including restraint system

271504001

Child EasyFix vacuum mattress Ready2go



### Accessories



GEAS260

Germa Vacuum Foot Pump