

## EasyFix vacuum mattress Ready2Go

Our brand new EasyFix Vacuum Mattress is made of flame-resistant PVC coated fabric. New inner bags with compartments for control of the granules, which means that it always has an even distribution of the content for a comfortable and stable fixation.

EASY FIX has a unique design that provides increased support of head, neck and back while packing volume has shrunk by 30% and weight by 20%. The top side of EasyFix has a reinforced base that covers the cervical spine. A unique wrapper function together with fixing straps gives a much better support than traditional systems. The development of this system is based on the X-spider functions where you think that the crossfixation provides clear benefits.

EasyFix is body shaped in design, providing a direct understanding of how to use the product while increasing the comfort and the ability to secure larger patients.

### Specifications:

Maximum shrinkage on length: 1%

Maximum shrinkage on width: 3%

Working temperatures 36 to + 75°C

### Material:

Outer material: Polyester fibre, PVC-coated (flame retardant)

Inner bag: Polyester

Granules: Styropor granules (flame retardant)

Treplate øvre del: Bjørk kryssfiner

### Other weigh information:

8,2 kg without fixation straps

8,5 kg with fixation straps

10 kg total weight including fixation straps + footrest pump

### Product information

- Head, neck and upper back part improved with wood part to secure stronger immobilisation
- New improved chambers system
- Light weighted and small storage
- Repair kit easy to use for effective repair
- Interchangeable valve
- Footrest high quality aluminium pump evacuate with more than 700 hpa
- Carry capacity of 220 kg, 10 handles
- Compliant to EN 1789 and 1865

### Productno.

271501000

EasyFix Vacuum Mattress

GE271501001

EasyFix vakuummattress Ready2Go



### Accessories



GEAS260

Germa Vacuum Foot  
Pump

