

## Ferno AS190 Orange Splint Set

The Ferno AS190 Vacuum splints are intended to stabilise injuries on patients during transport. They are suited for fixation of patients with hand, wrist, arm, leg, shoulder and/or hip injuries. The splints moulds around the individual fracture giving an optimal and comfortable support.

Set consists of: Arm, Long arm, Leg, Pump, Bag and Repair kit  
 AS100 and AS120 - Immobilisation of arm, forearm and wrist fractures.  
 AS140 - Immobilisation of lower leg on adult and full leg, hip or shoulder on child.  
 Also works as a vacuum mattress for babies.

Material:

Outer material: Polyester fibre, PVC coated (flame retardant)

Inner bag: Polyester

Granules: Styroporgranules (flame retardant)



### Product information

- Deliver consistent and rigid support to the patient
- Made of robust and efficient material to hold vacuum for more than 15 hours
- Multi fit valve
- Germa valve with "reverse air-lock" and supporting counter
- Hygienic design easy clean
- X-ray compatible and translucent
- Fits and can be attached to most stretchers on the market
- Repair kit easy to use for effective repair
- Interchangeable valve

### Productno.

GEAS120	Vacuum splint, long arm
GEAS140	Vacuum splint, long leg
GEAS100	Vacuum splint, short arm
GEAS190	Ferno Vacuum splint set in carrying bag



### Accessories



GEAS160

Germa Vacuum Hand  
Pump