

Ferno Orange Line Transfer Mattress

Safe Transfer Mattress has been developed in co-operation with our customers with focus on patient safety and comfort, this together with additional values of practical use to treat the patient.

Safe Transfer Mattress has been tested together with trauma teams to secure efficient handling plus quick and easy transfer of patients also high weight patients.

Safe Transfer Mattress has been crash tested and approved according to EN 1865 and EN 1789

Specification:

PVC-Polyester

4,5 cm polyeter + 1,5 cm Sensus (visko elastic). Total 60 mm

Handles: 12



Product information

- Equipped with unique rugged sliding material on the back of the mattress, enabling smooth transfer of mattress with patient from stretcher to stretcher, trolley or x-ray facility.
- Handles are reinforced to withstand the force required. This reinforcement also makes it possible to use the Germa Stretcher/Transfer Mattress in a patient lifting device.
- Ergonomic grip on the handles
- The mattress has a flexible inner plate to ensure correct fit on the stretcher.
- No matter how the stretcher is positioned, the mattress is flexible and will follow the curves of the stretcher.
- Robust welding - fully welded without any open seams or textiles, making hygienic cleaning of the mattress easy

Productno.

271022000

Stretcher/Transfer Mattress 16



Accessories



GE10001

Step by Step Evacuation Set