

Germa Baby Stroller Support

Germa Baby stroller support is intended to be placed inside a stroller for positioning of a child who have limited mobility and need support to sit comfortably and safely. With help of the vacuum technology the pillow is shaped exact after the child body and takes away every pleasure point and this make it possible to make it stiff via evacuation of air. This feature is to support the child in safe sitting positioning with good comfort for longer period.



Product information

- Shapes after the child's body
- Support in sitting position

Productno.

| | |
|-------------|-------------------------|
| 22218400000 | Germa Form Pillow 22x26 |
| 22218540000 | Germa Form Plus 22x34" |
| 22218340000 | Germa Form Plus 22x12" |
| 22218440000 | Germa Form Plus 22x26" |
| 22218600000 | Germa Form Pillow 16x20 |
| 22218140000 | Germa Form Plus 32x40" |
| 22218500000 | Germa Form Pillow 22x34 |
| 22218640000 | Germa Form Plus 16x20" |

Accessories



GEAS160

Germa Vacuum Hand Pump