

## Germa Blue Line Vacuum Pillows

These pillows are intended for customized, semi-permanent positioning anywhere contoured support is needed. Pillows are filled with small styrene granulates that reads your body shapes when you evacuate the air it will support strongly the chosen shape. When the valve is released, Germa Form Pillow regains its flexibility and is ready to be shaped again. Forms a firm, support for all parts of the body, reducing risk of pressure points. Use to accommodate structural deformities and help control posture. Shapes and reshapes in seconds. Easy to clean and disinfect.

Application of Germa Blue Vacuum Pillows provide customised support for a large range of different need and situations when long lasting fixation is desirable:

Semipermanent stabilisation  
Positioning  
Fixation

To support:

Stabilisation of children and adults with different types of disabilities  
Fixation at hand-, eye- and general surgery  
Tool for fixation at preliminary x-ray

When air is evacuated by use of e.i. the one-hand-operated vacuum pump, the pillow becomes stable in the chosen form and stiffness.  
When the valve is released, allowing air into the pillow, the Germa Vacuum Pillow regains its flexibility and is ready to be reshaped again.

Specification:

Germa Vacuum Pillows are made of PVC-film or PVC coated fabric with vacuum valve

Inner bag in nylon. Granules of polystyrene

Germa Vacuum Pumps are made of aluminium and robust plastic

The bigger foot rest vacuum pump has a higher degree of efficiency especially for larger vacuum pillows.



### Product information

- The Vacuum Pillows consist of an strong airtight shell containing larger quantity of granules made of polystyrene.
- Available in several sizes up to mattress size.
- Can be shaped and reshaped in seconds
- Density and firmness can be easily adjusted
- Durable, strong blue PVC-film in shell, easy to clean and disinfect for multiple client use
- 
- 

### Productno.

GE22396640000	40 x 50 cm Germa Vacuum Pillow, blue
GE22391140000	80 x 100 cm Germa Vacuum Pillow, blue

### Accessories



GEAS160

Germa Vacuum Hand  
Pump