

Germa Blue Splint Line

The Germa Vacuum splints are intended to stabilise injuries on patients during transport. They are suited for fixation of patients with hand, wrist, arm, leg, shoulder and/or hip injuries. The splints moulds around the individual fracture giving an optimal and comfortable support.

There are two different sets delivered:

1. Arm, Leg, Pump, Bag and Repair kit
2. Arm, Low Leg, Full Leg, Pump, Bag and Repair kit

Material:

Outer material: Polyester fibre, PVC coated (flame retardant)

Inner bag: Polyester

Granules: Styroporgranules (flame retardant)

Product information

- Working temperatures: - 36 to + 75°C
- Quality controlled production ISO 9001:2008
- 1 Year Warranty against production faults on production quality

Productno.



Accessories



GEAS160

Germa Vacuum Hand
Pump

