

Germa Military Splints MS MOD

The Germa Vacuum splints are intended to stabilise injuries on patients during transport. They are suited for fixation of patients with hand, wrist, arm, leg, shoulder and/or hip injuries. The splints moulds around the individual fracture giving an optimal and comfortable support.

Working temperatures 36 to + 75°C

Quality controlled production ISO 9001:2000

3 Year Warranty against production faults on production quality

Material

Outer material: Polyester fibre, PVC coated (flame retardant)

Inner bag: Polyester

Granules: Styroporgranules (flame retardant)

Set consists of: Arm, Low Leg, Full Leg, Pump, Bag and Repair kit

Product information

Productno.



Accessories



GEAS160 GEAS260
Germa Vacuum Pumps

