

Germa MOD Standard Line Mattress

The vacuum mattress is a combined stretcher and vacuum mattress. To save storage space the mattress can lie on the stretcher at all times, acting as a stretcher mattress. If vacuum immobilisation is needed, the vacuum mattress is integrated in the stretcher mattress and is ready for use immediately. Of course, the mattress can also be used as a stand-alone ordinary vacuum mattress.



Specifications

Maximum shrinkage on length: 1% Maximum shrinkage on width: 3%

Material

Outer material:

Polyester fibre, PVC coated (flame retardant)

Inner bag: Polyamide

Granules: Styropor granules (flame retardant)

Working temperatures -36 to +75°C

Product information

- Carry capacity 250 kg
- Approved according to European Standard EN 1865
- Quality controlled production ISO 9001:2000
- 3 Year Warranty against production faults on production quality

Productno.



Accessories



GEAS260

Germa Vacuum Foot Pump