

Germa Orange Line Splints

The Germa Orange Line is the original high quality long-life model.

The Germa Vacuum splints are intended to stabilise injuries on patients during transport. They are suited for fixation of patients with hand, wrist, arm, leg, shoulder and/or hip injuries. The splints moulds around the individual fracture giving an optimal and comfortable support.

Material:

Outer material: Polyester fibre, PVC coated (flame retardant)

Inner bag: Polyester

Granules: Styroporgranules (flame retardant)

Specifications:

Working temperatures 36 to + 75°C



Product information

- Immobilisation of arm, forearm and wrist fractures.
- Quality controlled production ISO 9001:2008
- 3 Year Warranty against production faults on production quality

Productno.



Accessories



GEAS160

Germa Vacuum Hand
Pump