

Mini Tranzfer Mattress

We have in cooperation with professional ambulance paramedics developed a very handy and useful tool for patient handling e.g. to move a patient from a stretcher to a bed or lift a patient from the floor or chair.

One of the big feature you will experience is when patient has to be rescued from a narrow environment e.g. bathroom. The great benefit is that you are able to bend the flexible stretcher/mattress around the patient and in step two lift and move the patient carrying in 10 possible strong handles.

To make usage as easy as possible, back of the Mini Trans mattress is made of friction reducing material to slide easy behind or under the patient. Glass fibre material on the back is very durable and can be used sliding on floor or ground.

Material:

Clear coloured green PVC-covered polyester

Lie area is with 10 mm thick Polyether foam



Product information

- Approved according to European Standard EN 1865.


Productno.

GE70050

Minitranz



Accessories



GE10000

Easy Lift carry support