

Mondial Mattress contoured

Mondial Mattress is developed in co-operation with Ferno International focusing on safety and comfort for the patient.

This mattress is designed with help of trauma teams to secure efficient handling, safety and comfort. Mattress is designed and construction is compliant to the European standard EN 1865 and EN 1789. This mattress comes with 2 year warranty, not covering for outer violence.



Product information

- Hygienic, all welded product, easy to clean with most common disinfectants
- Side support, contoured mattress to create safe feeling and comfort
- Designed to fit on top of Ferno Mondial stretchers, fixation with plastic hook. Important to use stretcher restraint system according to European standard EN 1789
- Velcro, hooks and loops adhesive opportunity for head part fixation

Productno.

GE23271664000

Mondial Mattress contoured

