

## Saver 190 Vacuum Splint Set

The Saver 190 Vacuum splint sets are intended to stabilise injuries on patients during transport. They are suited for fixation of patients with hand, wrist, arm, leg, shoulder and/or hip injuries. The splints moulds around the individual fracture giving an optimal and comfortable support.

Set consists of: Arm, Long arm, Leg, Pump, Bag and Repair kit

Working temperatures 36 to + 75°C

Quality controlled production ISO 9001:2008

1 Year Warranty against production faults on production quality

Material

Outer material: Polyester fibre, PVC coated (flame retardant)

Inner bag: Polyester

Granules: Styroporgranules (flame retardant)



### Product information

**Productno.**

005002032

Saver 100 Arm Splint

005002033

Saver 120 Arm Splint



### Accessories



GEAS160

Germa Vacuum Hand Pump