

Saver Blue Line Vacuum Mattress

We have developed the Blue Line / Saver mattress on the platform of our All In One success a vacuum mattress which is a combined stretcher and vacuum mattress. To save storage space the mattress can lie on the stretcher at all times, acting as a stretcher mattress. When vacuum immobilisation is needed, the vacuum mattress is integrated in the stretcher mattress and is ready for use immediately. Of course, the vacuum mattress can also be used as stand-alone ordinary vacuum mattress stored in the ambulance and covered by our bag including vacuum pump and repair kit.

Blue Line / SAVER is of simple inner chamber construction and do not have same control over the granulates and need extra preparation for good immobilisation before use.

Working temperatures 36 to + 75°C
Quality controlled production ISO 9001:2000
3 Year Warranty against production faults and production quality

Application:
The Saver is of high quality long-life model.

Specifications
Maximum shrinkage on length: 1%
Maximum shrinkage on width: 3%
Approved according to European Standard EN 1865

Material
Outer material: Polyester fiber, PVC-coated (flame retardant)
Inner bag: Polyamide
Granules: Styroporgranules (flame retardant)

NB! Delivered without straps. Straps are ordered separately!



Product information

- Vacuum mattress granulate-filled mattress without innerbag but with the possibility to refill the granules through the valve.
- Two different closing mechanisms on the straps, either Velcro or plastic buckle, allowing the user to choose his/her preference
- Valve with special designed hole, ensuring almost all suction devices fit directly into the valve, no adapters are needed. Made from plastic makes it ideal for X-ray.

Productno.



Accessories



GEAS250

Belt set 4 colours -
Velcro



GEAS260

Germa Vacuum Foot
Pump