

Slide/carry sheet 6 handles

Germa Carry Sheet contains a PVC fabric material welded with 6 carry handles on each side of the sheet.

This equipment is designed to support and help the paramedics to transfer a patient from a bed or floor on to a ambulance stretcher and in next step from the stretcher to a hospital bed. The purpose is to support patient with better comfort and paramedics will enjoy much appreciated ergonomic tool for the patient handling.

Material:

PVCpolyester fabric material (flame retardant)



Product information

- CEMarked
- X-ray and MRI translucent
- Approved according to European Standard EN 1865 and EN 1789

Productno.

GE24005000001

Slide/carry sheet 6 handles

