

Step by Step Funeral

The Step by Step has two belts with carabiners to connect to transfer mattress.

The system is used to slide a patients down stairs.

Step by Step is quick and easy to apply to the mattress and includes a restraint system for the patient to be secured on the mattress.
This solution deliver a secure evacuation of laying patients with ergonomics and less hassle for the rescue teams.

Product information

- Uses little space
- Weighs only 1,25 kg
- Can be used on all our transfer mattresses
- Can be used on weights up to 250 kg
- Easy to understand user manual

Productno.

GE15000080014

Step by Step Funeral



Accessories



GE271008800

Transfer Mattress Black