

## Vaccum Mattress All in One II P/V

Germa offers this second generation of All In One line of Vacuum mattresses: The Germa All in One II P/V is improved version from the original high quality long-life model All in One.

The vacuum mattress is a combined stretcher and vacuum mattress. To save storage space the mattress can lie on the stretcher at all times, acting as stretcher mattress. When vacuum immobilisation is needed, the vacuum mattress is integrated in the stretcher mattress and is ready for use immediately. Of course the mattress can also be used as a stand-alone ordinary vacuum mattress.

### Specifications

Maximum shrinkage on length: 1%

Maximum shrinkage on width: 3%

Working temperatures: -36 to +75°C

### Material

Outer material:

Polyester fibre, PVC-coated (flame retardant)

Inner bag: Polyamide

Granules: Styropor granules (flame retardant)

Wood board head part

Quality controlled production ISO 9001:2000

3 Year Warranty against production faults and production quality



## Product information

- Head and neck part improved with reinforcement for stronger immobilisation
- Safe use of 4-point restraint system from the stretcher according to EN 1789
- Improved chambers system to control granulates, quick and easy to prepare before use.
- Choice of three options of width 90, 100 or 115 cm capacity to give immobilisation for large patients
- This mattresses is crash tested by notified body for approval to the European standards.
- Compliant to EN 1789 and 1865

### Productno.

271303000

Vacuum Mattress All in One II P/V 90 cm

271301000

Vacuum Mattress All in One II P/V 115 cm

271302000

Vacuum Mattress All in One II P/V 100 cm



## Accessories



GEAS260

Germa Vacuum Foot Pump