

## Vacuummadrass All in one

All-in-One mattress has the same sturdy handle that transport the mattress and meets the EN1865 demands. The mattress can be used as transport mattress and vacuum the mattress as "one-piece".

All in One mattress is also excellent as a mobility aid such as the narrow stairs when the patient is both fixed and protected in a safe manner, and that all the robust handles give you that use multiple choice.

All-in-One mattress has the same sturdy handle that transport the mattress and meets the EN1865 demands. The mattress can be used as transport mattress and vacuum the mattress as "one-piece".

All in One mattress is also excellent as a mobility aid such as the narrow stairs when the patient is both fixed and protected in a safe manner, and that all the robust handles give you that use multiple choice.



### Product information

- Made of robust and efficient material to hold vacuum for more than 15 hours
- Extra layer built in under the mattress for safety and that the mattress last longer
- Carrying handles allow 2 people or more to carry which facilitates difficult handling
- Multi fit valve gives the possibility to use other available suction pumps
- Germa valve with "reverse air-lock" and supporting counter reduces hazzles.
- Safety belt fixation close to the mattress for safety and comfort
- Hygenic design maske it easier to clean
- Fits and can be attched to most stretchers on the market
- Repair kit which is easy to use for effective repair. Maximum shrinkage on length: 1%
- Maximum shrinkage on width: 3%
- Working temperatures 36 to + 75°C
- Approved according to European Standard EN 1865
- Quality controlled production ISO 9001:2000
- 3 Year Warranty against production faults and production quality

### Productno.

GE271203000	Germa All in one 90 cm wide, 16 handles
GE271201000	Germa All in one 115 cm wide, 16 handles
271104000	Germa All in One Line 80 cm wide, 12 handles
GE271103000	Germa All in One Line 90 cm wide, 12 handles
GE271102000	Germa All in one 100 cm wide, 12 handles
GE271101000	Germa All in One 115 cm wide, 12 handles
GE271202000	Germa All in one 100 cm wide, 16 handles



### Accessories



GEAS245

Bag for All in One II Vacuum mattress



GEAS250

Belt set 4 colours - Velcro



GEAS260

Germa Vacuum Foot Pump