



**AB GERMA®**

# Ferno F2 MD Comfort Stretcher/Cot Mattress

**Ferno F2 MB Mattress** is developed in co-operation with Ferno International focusing on safety and comfort for the patient on a Ferno F2 stretcher/cot.

This quality mattress is designed with help of trauma teams to secure efficient handling with safety and comfort for the patient. Mattress is designed to be compliant to the European standard EN 1865 and EN 1789.

### Features and benefits:

#### Safety

Easy to assembly and adjust stretcher restraint system according to standard regulation EN 1789

#### Comfort

Extra comfort by using a combination with Comfort foam which reads and support patients body and even out pressure to be able to deliver premium safe and comfort feeling when laying on the mattress.

#### Design

Designed to be adjusted in optimal position to bridge from stretcher/

- Hygienic, all welded product, easy to clean with most common disinfectants
- Patient support via Comfort foam in different density in the mattress to create safe feeling and comfort.
- Designed to fit on top of Ferno F2 stretcher/cot, fixation with hooks and loops.
- Important to use stretcher restraint system according to European standard EN 1789 and instructions.



Revers fixation with hooks and loops

Specification - Ferno F2 MD Stretcher /Cot Comfort Mattress	Dimensions and weight
Padova PVC- Polyester fibre flame retardant, compliant to EN 1865 all parts optimal density foam material	Length: 191 x 49 x 9 cm Weight: 4,1 kg
Ordering Specification	Art. No.
<b>Ferno F2 MD Stretcher/Cot Mattress</b> Ferno label	23271974000



**AB GERMA®**

**AB GERMA**  
P O Box 2127,  
Industrigatan 54-56  
SE-291 02 KRISTIANSTAD, SWEDEN  
Phone: +46-44 12 30 30