



**AB GERMA®**

## Heavy Duty Patient Handling

### Transfer Mattress Heavy Duty Patient Handling.

#### Application

Ferno Transfer Mattress Heavy Duty contains a PVC fabric material welded with 6 carry handles on each long side plus evacuation head and foot end handles of the mattress.

This equipment specifically designed to support and help the paramedics and rescuers to be able transfer and transport a heavy-duty patient. Designed to be able to slide and when needed carry a patient from an narrow environments and/or bed or floor on to a heavy-duty ambulance stretcher. Next step from the stretcher to a hospital bed. The purpose is to support handling of heavy-duty patient with safety and improved ergonomic, appreciated ergonomic tool for heavy-duty patient handling.

#### Features and Benefits

- 6 carry handles on each side make it easy to find balance for lifting when needed and carry a patient on board. Easy accessible handles up from ground lift.
- 2 handles on head and foot end to slide the patient on floors and staircases
- Sliding surface on the back of the Transfer Mattress
- Easy to clean fully welded product
- 2 Safety belts to secure the patient and cover sheet to protect and hold the patient
- Durable and strong material and design with a carry capacity of 380 kg
- Heavy Duty Transfer Mattress should be accessible when needed for challenging patient handling.



#### Specifications

Approved according to European Standard  
EN 1865  
CE—marked  
X-ray and MRI translucent

Length: 205 cm (80,70 in)  
Width: 143 cm (56,30 in)  
Weight: 11 kg (24,4 lb)

#### Material:

PVC—polyester fabric material (flame retardant)  
Rope: Nylon  
Foam: Sensus comfortable memory foam

**Ferno Tranfer Mattress Heavy-Duty**

**Art no: 271 050 000**

**203 cm long, 185cm wide, 60cm high  
Weight 10,5 kg**