



AB GERMA®

UNI EASY TRANSFER MATTRESS



PRODUCT SHEET

Uni Easy Transfer mattress has been developed in co-operation with professional paramedics with focusing on hygienic, ergonomic patient handling, patient safety and comfort.

Uni Easy Transfer mattress has been tested together with medical trained trauma teams to secure efficient handling plus quick and easy transfer of patients.

The Uni Easy Transfer mattress has been designed to fit most common stretchers available on the market.

UNI EASY TRANSFER MATTRESS



X-RAY MRI AND
CT TRANSLUCENT



FLAME RETARDANT



HF-VELDED SEAMS

PRODUCT SPECIFICATIONS

Approved according to European Standard EN 1865 and EN 1789
MDR 2017/745 Class 1 Medical device

PRODUCT INFO

- 5 handles on each side
- Glass fiber back side
- Carry capacity 250 Kg

Warranty period: 3 years

MATERIAL

- Outer material: Polyester fibre, PVC-coated, flame retardant.
- Inner material: 40 mm V - shaped Polyurethane Foam + 15 mm Sensus Foam

SPECIFICATION FOR THE UNI EASY TRANSFER MATTRESS

- Length: 190 cm
- Width: 50 cm
- Height: 6 cm
- Weight 4,7 kg

TEST SPECIFICATIONS

- Working temperature of the mattress is -36°C to +75°C

ORDER INFORMATION

Art no: 271030001
Uni Easy Transfer Mattress

Visiting address
AB Germa
Industrig. 54-56
SE-291 36 Kristianstad, Sweden

Post address
AB Germa
Box 2127
SE-291 02 Kristianstad, Sweden

www.germa.se
Phone: +46 (0)44 12 30 30
info@germa.se

BENEFITS

UNI EASY TRANSFER MATTRESS

Made to correspondent to EN 1789 and 1865

The Uni Easy Transfer mattress has been designed to fit most common stretchers available on the market.

High level of patient comfort due to the V-shaped foam plus the Sensus memory foam.

5 handles on each side, allow paramedics to easy lift, carry and slide the mattress with patient on it.

Sliding surface on the reverse of fiberglass material, reducing friction and allow sliding of the mattress e.g. on the floor/ground or from stretcher to stretcher or from stretcher to hospital bed.

For easier patient access during transport the handle part fold under the mattress on the stretcher.

It is very easy to apply the stretchers restraint system due to the easy access slits at the shoulder level on the Uni Easy Transfer mattress.

For patient safety the mattress should be secured on the stretcher with the stretchers restraint system.

It is hygienic as it is fully welded and easy to clean with most common disinfectants.