



AB GERMA®

# EASYFIX PLUS VACUUM SPLINTS



## PRODUCT SHEET

Germa EasyFix PLUS line of splints:

The EasyFix PLUS is built on the original high quality, long-life model.

EasyFix PLUS Vacuum splints are intended and designed to deliver premium stabilisation to injured patients during transport. They are suited for fixation of patients with hand, wrist, arm, leg, shoulder and/ or hip injuries. The splint moulds around the individual fracture giving an absolute optimal and very comfortable support to the patient.

EasyFix PLUS Vacuum splints has Anti-Bacterial hooks and loops as well as Improved hygienic BioSafe straps.

Designed and constructed to reduce needed preparation and is very reliable and a product you can trust, world class equipment.

# EASYFIX PLUS VACUUM SPLINTS



X-RAY MRI AND  
CT TRANSLUCENT



FLAME RETARDANT



MULTI FIT VALVE WITH LID



HF-WELDED SEAMS

## PRODUCT SPECIFICATIONS

Approved according to European Standard EN 1865 and EN 1789  
Compliant to MDR 2017/745 Class 1 Medical device

### PRODUCT INFO

- Anti-Bacterial hooks and loops
- Improved hygienic - BioSafe straps
- Flexible to positioning of the splint - premium immobilisation

Warranty period: 3 years

### MATERIAL

Outer material:

Polyester fibre, PVC coated (flame retardant)

Inner bag: Polyamide

Granules: Styropor granules (flame retardant)

### ORDER INFORMATION:

Art: Full Leg Splint

Art no: 272 111 000

Immobilisation of full leg, hip or shoulder injuries to adults.

Also works as a vacuum mattress for small children or babies.

Length: 136 cm

Width: 73 cm

Thickness: 3 cm

Weight: 1.55 kg



Art: Leg Splint

Art no: 272 112 000

Immobilisation of lower leg on adult and full leg, hip or shoulder on child. Also works as a vacuum mattress for babies.

Length: 98 cm

Width: 61 cm

Thickness: 3 cm

Weight: 1.15 kg



Art: Arm Splint

Art no: 272 113 000

Immobilisation of arm, forearm and wrist fractures.

Length: 74.5 cm

Width: 39 cm

Thickness: 3 cm

Weight: 0.75 kg



Art: Vacuum Pump

Art no: 272 000 701

One-hand-operated Vacuum Pump.

Weight: 0.6 kg



Visiting address

AB Germa

Industrig. 54-56

SE-291 36 Kristianstad, Sweden

Postal address

AB Germa

Box 2127

SE-291 02 Kristianstad, Sweden

[www.germa.se](http://www.germa.se)

Phone: +46 (0)44 12 30 30

[info@germa.se](mailto:info@germa.se)

## BENEFITS

Compliant to European standard EN 1865 and EN 1789  
MDR 2017/745 Class 1 Medical device

Unique inner bag with channel system, ensuring granulates are evenly dispersed. Quick and easy to use.

Robust weldings. Trustful.

PVC-coated fabric that is easy to clean and disinfect.

Full Leg splints with extra Velcro for strap allowing the splint to be used as vacuum mattress for babies and children.


Full leg splint is provided with additional Velcro straps to ensure one size fits all.

Velcro is sewn onto a piece of PVC coated fabric to reinforce the strength of the Velcro straps.

### ORDER INFORMATION SETS:

Art: Splint Set I  
Art no: 272 114 000  
Arm, Leg, Repair Kit, Pump  
and Bag  
Length: 75 cm  
Width: 37 cm  
Thickness: 18 cm  
Weight: 3.10 kg

Art: Splint Set II  
Art no: 272 115 000  
Arm, Leg, Full leg, Repair Kit, Pump  
and Bag  
Length: 75 cm  
Width: 37 cm  
Thickness: 18 cm  
Weight: 4.65 kg



Visiting address  
AB Germa  
Industrig. 54-56  
SE-291 36 Kristianstad, Sweden

Postal address  
AB Germa  
Box 2127  
SE-291 02 Kristianstad, Sweden

[www.germa.se](http://www.germa.se)  
Phone: +46 (0)44 12 30 30  
[info@germa.se](mailto:info@germa.se)

# EASYFIX PLUS VACUUM SPLINTS