



AB GERMA®

EASYFIX VACUUM SPLINTS



PRODUCT SHEET

Germa EasyFix PLUS line of splints:

The EasyFix PLUS is built on the original high quality, long-life model.

EasyFix PLUS Vacuum splints are intended and designed to deliver premium stabilisation to injured patients during transport. They are suited for fixation of patients with hand, wrist, arm, leg, shoulder and/ or hip injuries. The splint moulds around the individual fracture giving an absolute optimal and very comfortable support to the patient.

Designed and constructed to reduce needed preparation and is very reliable and a product you can trust, world class equipment.

EASYFIX VACUUM SPLINTS



X-RAY MRI AND
CT TRANSLUCENT



FLAME RETARDANT



MULTI FIT VALVE WITH LID



HF-WELDED SEAMS

PRODUCT SPECIFICATIONS

Approved according to European Standard
EN 1865 and EN 1789
Compliant to MDR 2017/745 Class 1 Medical device

PRODUCT INFO

- Flexible to positioning of the splint - premium immobilisation

Warranty period: 3 years

MATERIAL

Outer material:

Polyester fibre, PVC coated (flame retardant)

Inner bag: Polyamide

Granules: Styropor granules (flame retardant)

ORDER INFORMATION:

Art: Full Leg Splint

Art no: 272 101 000

Immobilisation of full leg, hip or shoulder injuries to adults.

Also works as a vacuum mattress for small children or babies.

Length: 136 cm

Width: 73 cm

Thickness: 3 cm

Weight: 1.55 kg



Art: Leg Splint

Art no: 272 102 000

Immobilisation of lower leg on adult and full leg, hip or shoulder on child. Also works as a vacuum mattress for babies.

Length: 98 cm

Width: 61 cm

Thickness: 3 cm

Weight: 1.15 kg



Art: Arm Splint

Art no: 272 103 000

Immobilisation of arm, forearm and wrist fractures.

Length: 74.5 cm

Width: 39 cm

Thickness: 3 cm

Weight: 0.75 kg



Art: Vacuum Pump

Art no: 272 000 701

One-hand-operated Vacuum Pump.

Weight: 0.6 kg



Visiting address

AB Germa

Industrig. 54-56

SE-291 36 Kristianstad, Sweden

Postal address

AB Germa

Box 2127

SE-291 02 Kristianstad, Sweden

www.germa.se

Phone: +46 (0)44 12 30 30

info@germa.se

BENEFITS

Compliant to European standard EN 1865 and EN 1789
MDR 2017/745 Class 1 Medical device

Unique inner bag with channel system, ensuring granulates are evenly dispersed. Quick and easy to use.

Robust weldings. Trustful.

PVC-coated fabric that is easy to clean and disinfect.

Full Leg splints with extra Velcro for strap allowing the splint to be used as vacuum mattress for babies and children.


Full leg splint is provided with additional Velcro straps to ensure one size fits all.

Velcro is sewn onto a piece of PVC coated fabric to reinforce the strength of the Velcro straps. This feature makes it far more durable.

ORDER INFORMATION SETS:

Art: Splint Set I
Art no: 272 104 000
Arm, Leg, Repair Kit, Pump
and Bag
Length: 75 cm
Width: 37 cm
Thickness: 18 cm
Weight: 3.10 kg

Art: Splint Set II
Art no: 272 105 000
Arm, Leg, Full leg, Repair Kit, Pump
and Bag
Length: 75 cm
Width: 37 cm
Thickness: 18 cm
Weight: 4.65 kg



Visiting address
AB Germa
Industrig. 54-56
SE-291 36 Kristianstad, Sweden

Postal address
AB Germa
Box 2127
SE-291 02 Kristianstad, Sweden

www.germa.se
Phone: +46 (0)44 12 30 30
info@germa.se

EASYFIX VACUUM SPLINTS



AB GERMA®

EASYFIX VACUUM SPLINTS



Full Leg Splint



Leg Splint



Arm Splint

