



AB GERMA®

**STRETCHER MATTRESS
MONDIAL WITH INTEGRATED
SIDESUPPORT**



PRODUCT SHEET

The Mondial Mattress is developed in co-operation with Ferno International focusing on safety and comfort for the patient.

This mattress is designed with help of trauma teams to secure efficient handling, safety and comfort.

The mattress has built in Sensus foam for extra comfort and stability of the patient.

STRETCHER MATTRESS MONDIAL WITH INTEGRATED SIDESUPPORT



X-RAY MRI AND
CT TRANSLUCENT



FLAME RETARDANT



HF-VELDED SEAMS

PRODUCT SPECIFICATIONS

Approved according to European Standard
EN 1865 and EN 1789
MDR 2017/745 Class 1 Medical device

PRODUCT INFO

Warranty period: 2 years

MATERIAL

- Outer material: Padova PVC- Polyester fibre. flame retardant
- Inner material: 80 mm V - shaped Polyurethane Foam + 15 mm Sensus Foam

SPECIFICATION FOR THE UNI EASY TRANSFER MATTRESS

- Length: 196 cm
- Width: 47 cm
- Height: 8 cm
- Weight 4,1 kg

TEST SPECIFICATIONS

- Working temperature of the mattress is -36°C to +75°C

ORDER INFORMATION

Art no: 271664001

Stretcher Mattress Mondial with integrated sidesupport

BENEFITS

Made to correspondent to EN 1789 and 1865
MDR 2017/745 Class 1 Medical devise

Easy to assembly and adjust stretcher
restraint system.

The PVC fabric is made strong and durable for
long and hard use.

The Mondial mattress has been designed to
fit the Ferno Mondial stretcher.

Side support via V-shaped foam in different
density in the mattress to create safe feeling
and comfort.

Extra pillow can be added to create added
comfort for patients head and neck part.

The pillow is designed to be adjusted in
position and to bridge when using head
extension part of the stretcher.

Fixation of the Mondial mattress on the
stretcher has two options, either with plastic
hook or with Velcro straps. PLEASE NOTE:
Important to always use stretcher restraint
system according to European standard EN
1789.

For patient safety the mattress should be
secured on the stretcher with the stretchers
restraint system.

It is hygenic as it is fully welded and easy to
clean with most common disinfectants.

STRETCHER MATTRESS MONDIAL WITH INTEGRATED SIDESUPPORT