



AB GERMA®

EASYFIX PLUS II VACUUM MATTRESS



Including our New 3 In 1 Carry Bag



PRODUCT SHEET

Germa innovative revolutionary designed line of Vacuum Mattresses:

EasyFix PLUS II is a development work based on earlier experiences from the original high quality materials traditionally used producing Germa world class vacuum mattresses.

EasyFix PLUS II vacuum mattress using improved body shape design to optimize the vacuum mattress capabilities. Offers 20% saving on storage space of the mattress compared to our All In One Vacuum Mattresses. EasyFix PLUS II unique construction which deliver improved fixation of upper back, neck when excluding cervical collar on the patient and reinforced base plate that cover all of cervical spine. This together with fixation restraint system using X-system to support and create the most effective immobilisation.

EASYFIX PLUS II VACUUM MATTRESS



REINFORCED WOOD
INNER SYSTEM



FLAME RETARDANT



HF-VELDED SEAMS



X-RAY MRI AND
CT TRANSLUCENT



MULTI FIT VALVE WITH LID



PRODUCT SPECIFICATIONS

Approved according to European Standard
EN 1865 and EN 1789
MDR 2017/745 Class 1 Medical device

PRODUCT INFO

- 5 handles on each side – carry capacity 250 Kg
- New improved chambers system
- Improved design of neck support to support patient without cervical collar
- Reinforced reverse of the mattress to reduce friction and improve possible handling
- Improve head strap system to have better grip
- Light weighted and small storage

Warranty period: 3 years

MATERIAL

- Outer material: Polyester fibre, PVC-coated,
- Flame retardant
- Inner bag: Polyester
- Granules: Polystyrene granules (flame retardant)
- Wood board upper part: Birch plywood

SPECIFICATIONS

- Loading capacity: 250 kg
- Length: 208 cm
- Width: 130/72 cm (Head part)
- Thickness: 4,5 cm
- Weight: 9,6 kg (excluding pump etc.)

TEST SPECIFICATIONS

- Maximum shrinkage on length: 1%
- Maximum shrinkage on width: 3%
- Working temperature of the mattress is -36°C to +75°C

ORDER INFORMATION

Art no: 271511100
Vacuum Mattress EasyFix PLUS II
Weight: 9,5 kg

Art no: 271511101
Vacuum Mattress EasyFix PLUS II
Ready2Go Kit
– including restraint system, footrest Vacuum Pump,
repair kit, bag
Weight: 12,5 kg

SPARE PARTS

Art no: 12633800000
Art: Germa vacuum valve Black

Art no: 15000000111
Art: EasyFix Repair Kit

Art no: 15000000034
Art: Easy Handle Set

Visiting address
AB Germa
Industrig. 54-56
SE-291 36 Kristianstad, Sweden

Post address
AB Germa
Box 2127
SE-291 02 Kristianstad, Sweden

www.germa.se
Phone: +46 (0)44 12 30 30
info@germa.se

EASYFIX PLUS II VACUUM MATTRESS

BENEFITS

Made to correspondent to EN 1789 and 1865
Approved for European standard and crash test for 10-G

Germa head pillow design including rigid wood part solution. Efficient, quick and improved neck and head immobilisation, with visible ears in one go.

Support upper back, neck and head for safe immobilization when not using cervical collar

One-piece back of the mattress, providing the patient with even spine support.

Design of the upper part allow the rescuer to place the arm(s) in two different options, make sure to give the rescuer access to patient arm for treatment requested via patient arm(s).

Reinforced reverse of the mattress to reduce friction.

Designed to allow improved patient handling.

Made from plastic makes it ideal for X-ray.
Allows preliminary X-ray diagnostics, direct in to x-ray no extra handling of patient

Germa improved chamber technology for quick and easy preparation.
Give rescuers quick handling and excellent immobilisation for the patient ensuring the granulate is evenly dispersed to give the patient stable and optimal immobilisation.

Handles with ergonomic grip, making patient lifts as comfortable as possible for rescuers.

Handles fixed to the bottom of mattress on piece, maintain rigid immobilisation when lifting the patient.

Robust welding—fully welded without any open seams or textiles, making hygienic cleaning of the mattress easy.

Valve with special designed hole, ensuring most suction devices fit directly into the valve, no adapters are needed.

Reverse air lock system, no closing of valve at stressful site.

Improved V-shaped design to make sure upper wide part will make all sized patients the very best immobilisation of the back, neck and head.

X-designed restraint belt fixation close to the patient, support immobilisation process and comfortable for patient. Supplied with green colour straps with fix lock, helping the rescuer to place the straps correctly and easy to stretch and adjust.

Outer material is made of flame retardant PVC-coated fabric.

PVC coated fabric makes it easy to clean and disinfect with most common disinfectants, easy to keep clean and ready for use.

Robust and efficient material to hold vacuum for more than 24 hours, safety for patient and make confiding immobilisation.

High load capacity 250 kg, safe to handle heavy duty patients.

Repair Kit easy to use for effective repair, quick to repair no return needed.

Folds easily, 20% less storage space required compared to most standard vacuum mattresses

High quality footrest aluminium vacuum pump evacuate with more than -700 hpa, durable and quick and time saving.

Visiting address
AB Germa
Industrig. 54-56
SE-291 36 Kristianstad, Sweden

Post address
AB Germa
Box 2127
SE-291 02 Kristianstad, Sweden

www.germa.se
Phone: +46 (0)44 12 30 30
info@germa.se