



AB GERMA®

VACUUM SPLINTS AS LINE



PRODUCT SHEET

The AS 190 Line of splints:

The AS 190 Line is Germa's Low-price range of splints, however, regardless of lower price they still have Germa's high quality and long-life.

The AS 190 Line splints are intended and designed to deliver stabilisation to injured patients during transport. They are suited for fixation of patients with hand, wrist, arm, leg, shoulder and/or hip injuries.

The splint moulds around the individual fracture giving an absolute optimal and very comfortable support to the patient.

VACUUM SPLINTS AS LINE



X-RAY MRI AND
CT TRANSLUCENT



FLAME RETARDANT



MULTI FIT VALVE WITH LID



HF-WELDED SEAMS

PRODUCT SPECIFICATIONS

Approved according to European Standard
EN 1865 and EN 1789
Compliant to MDR 2017/745 Class 1 Medical device

PRODUCT INFO

Warranty period: 3 years

MATERIAL

Outer material:
Polyester fibre, PVC coated (flame retardant)
Granules: Styropor granules (flame retardant)

ORDER INFORMATION:

Art: Vacuum Splint AS 140
Art no: 005002014
Immobilisation of full leg, hip
or shoulder injuries to adults.
Also works as a vacuum mattress
for small children or babies.
Length: 112 cm
Width: 70 cm
Thickness: 3 cm
Weight: 1.15 kg



Art: Vacuum Splint AS 120
Art no: 005002013
Immobilisation of lower leg
on adult and full leg, hip or
shoulder on child. Also works as
a vacuum mattress for babies.
Length: 75 cm
Width: 63 cm
Thickness: 3 cm
Weight: 0,75 kg



Art: Vacuum Splint AS 100
Art no: 005002012
Immobilisation of arm,
forearm and wrist fractures.
Length: 55 cm
Width: 36 cm
Thickness: 3 cm
Weight: 0.75 kg



Art: Vacuum HandPump
Art no: 272000701
One-hand-operated
Vacuum Pump.
Weight: 0.6 kg



Visiting address
AB Germa
Industrig. 54-56
SE-291 36 Kristianstad, Sweden

Postal address
AB Germa
Box 2127
SE-291 02 Kristianstad, Sweden

www.germa.se
Phone: +46 (0)44 12 30 30
info@germa.se

BENEFITS

Compliant to European standard EN 1865 and EN 1789
MDR 2017/745 Class 1 Medical device

Robust weldings. Trustful.

PVC-coated fabric that is easy to clean and disinfect.

Full Leg splints with extra Velcro for strap allowing the splint to be used as vacuum mattress for babies and children.

Germa valve with "reverse air-lock".

X-ray compatible and translucent

VACUUM SPLINTS AS LINE

ORDER INFORMATION SETS:

Art: Vacuum Splint Kit AS190
Art no: 005002210
Short arm, long arm, long Leg,
Repair Kit, Pump and Bag
Length: 82 cm
Width: 32 cm
Thickness: 18 cm
Weight: 4.30 kg



Visiting address
AB Germa
Industrig. 54-56
SE-291 36 Kristianstad, Sweden

Postal address
AB Germa
Box 2127
SE-291 02 Kristianstad, Sweden

www.germa.se
Phone: +46 (0)44 12 30 30
info@germa.se