



AB GERMA®

VACUUM SPLINTS BLUE LINE



PRODUCT SHEET

Germa Blue Line of splints:

The Blue Line is Germa's Low-price range of splints, however, regardless of lower price they still have Germa's high quality and long-life.

The Blue Line splints are intended and designed to deliver stabilisation to injured patients during transport. They are suited for fixation of patients with hand, wrist, arm, leg, shoulder and/or hip injuries.

The splint moulds around the individual fracture giving an absolute optimal and very comfortable support to the patient.

VACUUM SPLINTS BLUE LINE



X-RAY MRI AND
CT TRANSLUCENT



FLAME RETARDANT



MULTI FIT VALVE WITH LID



HF-WELDED SEAMS

PRODUCT SPECIFICATIONS

Approved according to European Standard
EN 1865 and EN 1789
Compliant to MDR 2017/745 Class 1 Medical device

PRODUCT INFO

Warranty period: 3 years

MATERIAL

Outer material:
Polyester fibre, PVC coated (flame retardant)
Granules: Styropor granules (flame retardant)

ORDER INFORMATION:

Art: Vacuum Splint Full Leg
Art no: 272701000
Immobilisation of full leg, hip
or shoulder injuries to adults.
Also works as a vacuum mattress
for small children or babies.
Length: 136 cm
Width: 73 cm
Thickness: 3 cm
Weight: 1.55 kg



Art: Vacuum Splint Lower leg
Art no: 272702000
Immobilisation of lower leg
on adult and full leg, hip or
shoulder on child. Also works as
a vacuum mattress for babies.
Length: 98 cm
Width: 61 cm
Thickness: 3 cm
Weight: 1.15 kg



Art: Vacuum Splint Forearm
Art no: 272707000
Immobilisation of arm,
forearm and wrist fractures.
Length: 74.5 cm
Width: 39 cm
Thickness: 3 cm
Weight: 0.75 kg



Art: Vacuum HandPump
Art no: 272000701
One-hand-operated
Vacuum Pump.
Weight: 0.6 kg



Visiting address
AB Germa
Industrig. 54-56
SE-291 36 Kristianstad, Sweden

Postal address
AB Germa
Box 2127
SE-291 02 Kristianstad, Sweden

www.germa.se
Phone: +46 (0)44 12 30 30
info@germa.se

BENEFITS

Compliant to European standard EN 1865 and EN 1789
MDR 2017/745 Class 1 Medical device

Robust weldings. Trustful.

PVC-coated fabric that is easy to clean and disinfect.

Full Leg splints with extra Velcro for strap allowing the splint to be used as vacuum mattress for babies and children.

VACUUM SPLINTS BLUE LINE

ORDER INFORMATION SETS:

Art: Vacuum Splint Set Blue Line
Art no: 272704000
Arm, Leg, Repair Kit, Pump and Bag
Length: 75 cm
Width: 37 cm
Thickness: 18 cm
Weight: 3.10 kg

Art: Vacuum Splint Set Blue Line +
Full leg Splint
Art no: 272705000
Arm, Leg, Full leg, Repair Kit, Pump
and Bag
Length: 75 cm
Width: 37 cm
Thickness: 18 cm
Weight: 4.65 kg



Visiting address
AB Germa
Industrig. 54-56
SE-291 36 Kristianstad, Sweden

Postal address
AB Germa
Box 2127
SE-291 02 Kristianstad, Sweden

www.germa.se
Phone: +46 (0)44 12 30 30
info@germa.se