



AB GERMA®

- ◇ MINI TRANSFER MATTRESS
- ◇ MINI TRANSFER MATTRESS MOD



PRODUCT SHEET

Together with professional paramedics and ergonomic specialist we developed a handy and useful mini transfer for patient handling. The objective was to reduce the hassle when bringing up a patient from the floor, bath tub, chair or move a patient from a bed over to a stretcher.

One big benefit is that you are able to bend the stretcher/mattress around or under the patient, in step two slide or lift the patient make use of the 10 available strong ergonomic handles.

Mini Transfer is suitable to be used in combination with Basket Stretcher to ease of patient handling in and out from the stretcher. In addition Mini Transfer deliver comfort to the patient.

- ◇ MINI TRANSFER MATTRESS
- ◇ MINI TRANSFER MATTRESS MOD



X-RAY MRI AND
CT TRANSLUCENT



FLAME RETARDANT



HF-VELDED SEAMS

PRODUCT SPECIFICATIONS

MDR 2017/745 Class 1 Medical device

PRODUCT INFO

- 4 handles on each side and 2 handles at the head.
- Foam filled patient area
- Glass fiber back side to protect and slide easier
- Max load capacity: 250 kg

Warranty period: 3 years

MATERIAL

- Outer material: Polyester fibre, PVC-coated, flame retardant
- Inner material: 10 mm thick Polyurethane Foam

SPECIFICATION FOR THE SPLINT SET

- Length: 136 cm
- Width: 78 cm
- Height: 1,5 cm
- Weight 1,8 kg

TEST SPECIFICATIONS

- Maximum shrinkage on length: 1%
- Maximum shrinkage on width: 3%
- Working temperature of the mattress is -36°C to +75°C

ORDER INFORMATION

Art no: 271040000

MINI TRANSFER MATTRESS

ORDER INFORMATION

Art no: 271048000

MINI TRANSFER MATTRESS MOD

SPARE PARTS

Art no: 15000000121

Art: REPAIR KIT MATTRESSES BLACK/GRAY

SPARE PARTS

Art no: 15000000124

Art: REPAIR KIT TRANSFER MOD

Visiting address
AB Germa
Industrig. 54-56
SE-291 36 Kristianstad, Sweden

Post address
AB Germa
Box 2127
SE-291 02 Kristianstad, Sweden

www.germa.se
Phone: +46 (0)46 12 30 30
info@germa.se

BENEFITS

Made to correspondent to EN 1789 and 1865

One big feature you will experience using Mini Transfer is that it fits great in a basket stretcher, the patient can then lay more comfortable and safer in the basket stretcher.

When patient has to be rescued from a narrow or difficult environment e.g. the forrest, the smaller size makes it easy to use and to manouver around.

The positions of the handles makes it possible to carry/move the patient with great flexibility.

To make usage easy as possible, back of the Mini Trans mattress is made with friction reducing material to slide easy behind or under the patient.

Glass fibre material on the back is very durable and can be used sliding on floor or ground.

- ◇ MINI TRANSFER MATTRESS
- ◇ MINI TRANSFER MATTRESS MOD