



**AB GERMA®**

**TRANSFER MATTRESS  
BLUE LINE 12 HANDLES**



# PRODUCT SHEET

Germa's Transfer Mattress Blue Line is intended to be used both as a comfort mattress for the animal when laying on the examination table, and as a transfer mattress when the animal has to be transferred from the examination table to, for example an X-Ray area. You don't need to lift the animal, you lift and move with the transfer mattress. It is similar to Germa's Veterinary Transfer Mattress but bigger so it can accommodate larger animals.

# TRANSFER MATTRESS BLUE LINE 12 HANDLES



X-RAY MRI AND  
CT TRANSLUCENT



FLAME RETARDANT



HF-VELDED SEAMS

## PRODUCT SPECIFICATIONS

### PRODUCT INFO

- 6 handles on each side
- Glass fiber back side
- Carry capacity 250 Kg

Warranty period: 3 years

### MATERIAL

- Outer material: Polyester fibre, PVC-coated, flame retardant.
- Inner material: 60 mm thick Polyurethane Foam

### SPECIFICATION FOR THE SPLINT SET

- Length: 195 cm
- Width: 50 cm
- Height: 6,5 cm
- Weight 5,3 kg

### TEST SPECIFICATIONS

- Maximum shrinkage on length: 1%
- Maximum shrinkage on width: 3%
- Working temperature of the mattress is -36°C to +75°C

### ORDER INFORMATION

Art no: 271708000

Transfer Mattress Blue Line 12 Handles

## BENEFITS

6 handles on each side, allow many people to lift and carry.

Handle part fold under the mattress.

Sliding surface on the reverse of fiberglass material, reducing friction and allow sliding of the mattress e.g. from examination table to another examination table.

Hygienic full welded easy to clean with most common disinfectants.

Visiting address  
AB Germa  
Industrig. 54-56  
SE-291 36 Kristianstad, Sweden

Post address  
AB Germa  
Box 2127  
SE-291 02 Kristianstad, Sweden

www.germa.se  
Phone: +46 (0)44 12 30 30  
info@germa.se