



AB GERMA®

VETERINARY TRANSFER MATTRESS



PRODUCT SHEET

Germa's Veterinary Transfer Mattress is intended to be used both as a comfort mattress for the animal when laying on the examination table, and as a transfer mattress when the animal has to be transferred from the examination table to, for example an X-Ray area. You don't need to lift the animal, you lift and move with the transfer mattress.

VETERINARY TRANSFER MATTRESS



X-RAY MRI AND
CT TRANSLUCENT



FLAME RETARDANT



HF-VELDED SEAMS

PRODUCT SPECIFICATIONS

PRODUCT INFO

- 3 handles on each side
- Glass fiber back side
- Carry capacity 180 Kg

Warranty period: 3 years

MATERIAL

- Outer material: Polyester fibre, PVC-coated, flame retardant.
- Inner material: 60 mm Polyurethane Foam

Specification for the Uni Easy Transfer Mattress

- Length: 130 cm
- Width: 50 cm
- Height: 6 cm
- Weight 3,7 kg

TEST SPECIFICATIONS

- Working temperature of the mattress is -36°C to +75°C

ORDER INFORMATION

Art no: 21271050000

Transfer Mattress Veterinarian

Visiting address
AB Germa
Industrig. 54-56
SE-291 36 Kristianstad, Sweden

Post address
AB Germa
Box 2127
SE-291 02 Kristianstad, Sweden

www.germa.se
Phone: +46 (0)44 12 30 30
info@germa.se