



**AB GERMA®**

# WHEEL CHAIR & STROLLER SUPPORT PILLOWS



# WHEEL CHAIR & STROLLER SUPPORT PILLOWS

## PRODUCT SHEET

Germa Form wheel chair support is intended to be placed inside a wheel chair or a stroller or for that matter in rehabilitation pools to help positioning a child who have limited mobility and need support to sit comfortably and safely.

With help of the vacuum technology the pillow is shaped after the child's body to support to sit in a good position. The pillow also takes away possible pressure point and makes it more comfortable for the child.

## PRODUCT SPECIFICATIONS

### PRODUCT INFO

- Is used to accommodate structural deformities and help control posture.
- Can be used in water for support in rehabilitation bathtubs and whirlpools.
- Easy to clean and disinfect.
- All models are delivered with repair kit and manual

### MATERIAL

- Outer material: Polyester fibre, PVC-coated, flame retardant
- Inner bag: Nylone
- Granules: Polystyrene

### SPECIFICATIONS

- Varies depending on size of the Wheel Chair Pillow

### ORDER INFORMATION

Art no: 22219040004

Positioning Pillow Wheelchair XL

Art no: 22219040003

Positioning Pillow Wheelchair L

Art no: 22219040002

Positioning Pillow Wheelchair M

Art no: 22219040001

Positioning Pillow Wheelchair S

---

#### Visiting address

AB Germa  
Industrig. 54-56  
SE-291 36 Kristianstad, Sweden

#### Post address

AB Germa  
Box 2127  
SE-291 02 Kristianstad, Sweden

[www.germa.se](http://www.germa.se)

Phone: +46 (0)44 12 30 30  
[info@germa.se](mailto:info@germa.se)