



AB GERMA®

STRETCHER MATTRESS
VIPER WITH
SIDESUPPORT



PRODUCT SHEET

The Viper Mattress is developed in co-operation with Ferno International focusing on safety and comfort for the patient.

This mattress is designed with help of trauma teams to secure efficient handling, safety and comfort.

The mattress has sidesupport for extra comfort and stability of the patient.

STRETCHER MATTRESS VIPER WITH SIDESUPPORT



X-RAY MRI AND
CT TRANSLUCENT



FLAME RETARDANT



HF-VELDED SEAMS

PRODUCT SPECIFICATIONS

Approved according to European Standard
EN 1865 and EN 1789
MDR 2017/745 Class 1 Medical device

PRODUCT INFO

Warranty period: 2 years

MATERIAL

- Outer material: Padova PVC- Polyester fibre. flame retardant
- Inner material: 40/60/100 mm Polyurethane Foam + 15 mm Sensus Foam

SPECIFICATION FOR THE UNI EASY TRANSFER MATTRESS

- Length: 199 cm
- Width: 47 cm
- Height: 4/6/10 cm
- Weight 4,0 kg

TEST SPECIFICATIONS

- Working temperature of the mattress is
-36°C to +75°C

ORDER INFORMATION

Art no: 271663000
Stretcher Mattress Viper

Visiting address
AB Germa
Industrig. 54-56
SE-291 36 Kristianstad, Sweden

Post address
AB Germa
Box 2127
SE-291 02 Kristianstad, Sweden

www.germa.se
Phone: +46 (0)44 12 30 30
info@germa.se

BENEFITS

Made to correspondent to EN 1789 and 1865
MDR 2017/745 Class 1 Medical devise

Easy to assembly and adjust stretcher
restraint system.

The PVC fabric is made strong and durable for
long and hard use.

The Viper mattress has been designed to fit
the Ferno stretcher.

Extra comfort by using a combination with
Sensus foam which reads and support
patients body and even out pressure to
deliver premium safe and comfort feeling
when laying on the mattress.

Extra pillow can be added to create added
comfort for patients head and neck part.

The pillow is designed to be adjusted in
position and to bridge when using head
extension part of the stretcher.

Fixation of the Viper mattress on the stretcher
is made by hook and Loops. PLEASE NOTE:
Important to always use stretcher restraint
system according to European standard EN
1789.

For patient safety the mattress should be
secured on the stretcher with the stretchers
restraint system.

It is hygenic as it is fully welded and easy to
clean with most common disinfectants.

STRETCHER MATTRESS VIPER WITH SIDESUPPORT