




AB GERMA®

**TRANSFER MATTRESS
12 HANDLES**



PRODUCT SHEET

Transfer Mattress 12 Handles has been used in the ambulance service for many years with very good experience of reducing the hassle, decrease forces required in evacuation and patient handling, ergonomic and handy. The construction is very durable and strong so it can be used in ruff terrain such as a forrest and other outdoor areas and not only on smooth indoor areas.

The Transfer Mattress also has a beneficial use as a transfer mattress on a stretcher to give comfort/decrease pressure wounds. The mattress is thinner than a regular stretcher mattress but still gives enhanced comfort for the patient.

TRANSFER MATTRESS 12 HANDLES



X-RAY MRI AND
CT TRANSLUCENT



FLAME RETARDANT



HF-VELDED SEAMS

PRODUCT SPECIFICATIONS

Approved according to European Standard EN 1865 and EN 1789
MDR 2017/745 Class 1 Medical devise

PRODUCT INFO

- 6 handles on each side
- Glass fiber back side

Warranty period: 3 years

MATERIAL

- Outer material: Polyester fibre, PVC-coated, flame retardant.
- Inner material: 30 mm thick Polyurethane Foam

SPECIFICATION

- Length: 187 cm
- Width: 65 cm
- Height: 3,0 cm
- Weight 4,5 kg

ORDER INFORMATION

Art no: 271008800
TRANSFER MATTRESS 12 HANDLES Black

SPARE PARTS

Art no: 15000000121
Art: REPAIR KIT MATTRESSES BLACK

Visiting address
AB Germa
Industrig. 54-56
SE-291 36 Kristianstad, Sweden

Post address
AB Germa
Box 2127
SE-291 02 Kristianstad, Sweden

www.germa.se
Phone: +46 (0)44 12 30 30
info@germa.se

BENEFITS

TRANSFER MATTRESS 12 HANDLES

Made to correspondent to EN 1789 and 1865

Easy to carry or slide over a body on the mattress with 12 strong handles, 6 on each long side.

Handles part fold under the mattress on the stretcher.

Sliding surface on the reverse of fiberglass material, reducing friction and allow sliding of the mattress e.g. on the floor/ground or from stretcher to stretcher.

Should be used with restraints to secure the body on top for safety reasons.

Comfort, ethical and ergonomic handling of the body.

Hygienic full welded easy to clean with most common disinfectants.