



AB GERMA®

GERMA VACUUM SPLINT SET MOD



PRODUCT SHEET

Germa line of Vacuum Splints MOD:

The Germa Vacuum Splints MOD's are built on the original high quality, long-life models that you are used to.

Germa Vacuum Splints MOD are intended and designed to deliver premium stabilisation to injured patients during transport. They are suited for fixation of patients with hand, wrist, arm, leg, shoulder and/or hip injuries.

The splint moulds around the individual fracture giving an absolute optimal and very comfortable support to the patient.

Designed and constructed to reduce needed preparation and is very reliable in it's function.

All Germa MOD products are made in a Non IR reflective PVC.

GERMA VACUUM SPLINT SET MOD



PRODUCT SPECIFICATIONS

The product is in conformity with Council MDR 2017/745 Class 1 concerning Medical Devices.

PRODUCT INFO

- Wrist Splint
- Arm Splint
- leg Splint
- Full Leg Splint
- Repair kit
- Hand operated Pump
- Carry bag.

Warranty period: 3 years

MATERIAL

- Outer material: Polyester fibre, PVC-coated,
- flame retardant and Non IR reflective
- Inner bag: Polyester
- Granules: Polystyrene granules (flame retardant)
- Carry Bag: Polyester

SPECIFICATION FOR THE SPLINT SET

- Length: 83 cm
- Width: 32 cm
- Height: 32 cm
- Weight 4,75 kg

TEST SPECIFICATIONS

- Maximum shrinkage on length: 1%
- Maximum shrinkage on width: 3%
- Working temperature of the mattress is -36°C to +75°C

ORDER INFORMATION

Art no: 272005800

GERMA VACUUM SPLINT SET+FULL LEG VACUUM
SPLINT MOD

– including hand operated Vacuum Pump, repair kit,
bag

SPARE PARTS

Art no: 15005002817

Art: Repair Kit Olive Green

Art no: 12633800000

Art: Germa vacuum valve Black

Visiting address
AB Germa
Industrig. 54-56
SE-291 36 Kristianstad, Sweden

Post address
AB Germa
Box 2127
SE-291 02 Kristianstad, Sweden

www.germa.se
Phone: +46 (0)44 12 30 30
info@germa.se

BENEFITS

Made to correspondent to EN 1789 and 1865

GERMA VACUUM SPLINT SET MOD

Unique inner bag with channel system, ensuring granules are evenly dispersed and quicker to use.

Robust welding, ensuring long life quality.

Leg splints with extra hooks and loops for strap allowing the splint to be used as vacuum mattress for babies and children.

Full leg splint is provided with additional hooks and loops straps to ensure one size fits all.

Hooks and loops is sewn onto a piece of PVC coated fabric to reinforce the strength of the Velcro straps.

High volume pump which is quick to use. Efficient so it saves time.

Multi fit valve.

Use almost any available pump.