




AB GERMA®

- ◇ TRANSFER SHEET UK
- ◇ BARIATRIC TRANSFER SHEET



PRODUCT SHEET

This equipment is designed to support and help the paramedics to transfer a patient from a bed or floor on to a ambulance stretcher and in next step from the stretcher to a hospital bed. The purpose is to support patient handling for the paramedics who will enjoy much appreciated ergonomic tool for patient handling and evacuation from narrow environment where you cannot go with a stretcher or cot.

◇ TRANSFER SHEET UK ◇ BARIATRIC TRANSFER SHEET



FLAME RETARDANT



X-RAY MRI AND
CT TRANSLUCENT



HF-VELDED SEAMS

PRODUCT SPECIFICATIONS

Approved according to EU Standard EN 1865.
CE Marked.

PRODUCT INFO

- 6 handles on each side make it easy to find balance for lifting and carry a patient.
- The bariatric versions handles are re-inforced to handle higher weight.

Warranty period: 2 years

MATERIAL

- PVC: Polyester fabric material (flame retardant)
- Rope: Nylon (in Bariatric only)

TRANSFER SHEET UK

SPECIFICATIONS

- Length: 185 cm
- Width: 80 cm (97 cm incl. handles)
- Weight 1,5 kg

- Carry capacity: 200 kg

ORDER INFORMATION

Art no: 005000004

Transfer Sheet UK

BARIATRIC TRANSFER SHEET

SPECIFICATIONS

- Length: 192 cm
- Width: 95,5 cm (113 cm incl. handles)
- Weight 1,9 kg

- Carry capacity: 380 kg

ORDER INFORMATION

Art no: 005000008

Bariatric Transfer Sheet

BENEFITS

- 6 handles on each side make it easy to find balance for lifting and carry a patient on board
- Easy to clean fully welded product.
- Durable and strong material and design
- X-ray, CT and MRI translucent.

Visiting address
AB Germa
Industrig. 54-56
SE-291 36 Kristianstad, Sweden

Post address
AB Germa
Box 2127
SE-291 02 Kristianstad, Sweden

www.germa.se
Phone: +46 (0)44 12 30 30
info@germa.se